Item Database File Download and Import Instructions:

Step 1: Click on the download link

Step 2: Select the Run option



Step 3: Select the **Run** option. Ignore the message regarding Unknown Publisher. This is referring to the WinZip Self-Extractor utility and not the Favorites Folders.



Step 4: Select the **OK** option to automatically extract the file(s) into the VeraDyne\Avatar 4.0\Import Folder

Physica Energetix Test Kits Import	X	3
Select OK to automatically Self-Extrac Files\VeraDyne\VeraDyne Avatar 4.0\J Remember you will need to use the " Menu-bar option) to Import this file.	t this file to the C:\Program import folder. Import Items" option (found in the Items	
	OK Cancel	

Step 5: Open Avatar 4.0 and select the Item Test Toolbar-button

<u>Session Edit</u> Items <u>Points</u> <u>H</u> oldTrays <u>TraceHist</u>	ory
📙 📴 Item Test 🛛 👋 Point Test 🛛 🔌 Trace	•
Active Trace Quick View Item Viewer	
100-0.00	
0.00	
90-	

Step 6: Select Items from the Menu-bar and then select Import Items



Step 7: The software will automatically open the contents of the Import folder displaying the Item Database file that you have downloaded. From the list select the folder that you wish to import then select **Open**

Select source	database:				X
Look <u>i</u> n:	import		• +	- 🖬 📩	
Ca.	Name	*		Date modified	Туре
Recent Places	Physica En	ergetix Test Kits.mdb		3/24/2009 5:27 PM	Micr
Desktop					
Robert Eanes					
Computer					
<u>_</u>					
Network					
	•	III			P.
	File <u>n</u> ame:	Physica Energetix Test Kit	s.mdb	▼ Op	en
	Files of type:	Database Files (*.mdb)		▼ Car	ncel

Step 8: The Avatar Item Import dialog-box will open displaying the contents of the database file (the list of items in the file). Select the **Import All Items** option-button

iource database: :\Program Files\VeraDyne\Ve	eraDyne Avatar 4.0\impoi	t\Physica Energetix Tesl	t Kits.mdb	🔍 Browse	
otal items: 653 Selected items	ε O		_		
Item Name			▲	mport All Items	
60 Hz AC Emission					
A/Fujian/411/2002			Imp	ort Selected Items	
A/Wyoming/03/2003 (H3N2)					
Acid Etch				Select All	
Acinic Cell Carcinoma					
Acquired Immunity				Select None	
Acrylic resins					
ACTH				Invert Selection	
Acute Herpetic Gingivostoma	titis				
Acute Necrotising Ulcerative	Icerative Gingivitis		S	elect Attributes	
Acute Necrotizing Ulcerative	Gingivitis				
Acute Osteomyelitis					
Acute Pulpitis			•		
•		•			
Clear existing attribute-state	s when undating duplica	te items			
elear entering attribute elate	io mich apadang dapilod	e nome.			
				Close	

A message-box window will display the progress of the import

Database Operation
Importing 653 items
174: Dendritic Cell
Cancel

Step 9: A message will appear verifying a successful import. Select **OK** to finish the process.



After completing the import, close the Avatar Item Import dialog-box. You will need to reload the entire Item database before the import will take effect. To do this you will need to close the Avatar software completely. Skip the backup step when you exit. Wait about 20 seconds and then open the Avatar software again. Remember to select the **Item Test** button and you should be able to see your newly import database items in the **Item Tree**.